



# VIVID Super PS Developer

Positive Plate Developer for PS plates

## Product Usage

---

Developer Concentrate for Positive Analogue Plates - ready to use  
universal for all popular positive working PS plates

## Characteristics

---

- Highly concentrated (1 part PS Developer + 9 parts water) - cost efficient
- Good buffer capacity in the plate processor - long lasting
- No or very low residues in the processor tank
- Universal in use - for all popular positive working plates

## Application

---

Dilute 1 part PS Developer Concentrate with 9 parts of water and mix well.

for replenishing dilute 1 part Super PS Concentrate with 6 parts of water

processor adjustments: Temperature: 22°C-26°C, start with 24°C

Processor speed: 20-40 seconds, start with 30 seconds

Replenishing rate: 15-40ml/m<sup>2</sup>, start with 20ml/m<sup>2</sup>

## Packaging

---

5/20 liter

---

**Vivid Manufacturing Co. Pvt. Ltd., 14/7, Site-IV, Industrial Area, Sahibabad-201010,U.P.**  
E-mail : [info@vividchem.com](mailto:info@vividchem.com) WebSite: [www.vividchem.com](http://www.vividchem.com)